



The next generation instruments in Hysteroscopy

The new and versatile solution for office and surgical Hysteroscopy. VERSATIP combines all the surgeon's needs with the requirements for appropriate cleaning in the CSSD.

- ▶ Detachable – one handle for all inserts
- ▶ Color code – the easy identification of the tip
- ▶ Two handles – self holding and standard
- ▶ The easy connection – the reliable ball tip technique
- ▶ Light weight handle – the perfect fit for your fingers
- ▶ Make cleaning easy – flushing port for the complete instrument
- ▶ 360° rotation - just with your finger tip
- ▶ All parts can be ordered separately



PATENTED
10 2022 107 371.8






Multiple tips design



Color coded inserts



Unique, detachable system



Flushing port
for easy cleaning




Two handles -
multiple options

FORCEPS & PUNCHES


- ▶ Automatic self-holding
- ▶ Black color coded



351-100-550
Self-holding handle for
forceps and punches
(recommended handle)



Grasping Forceps with toothed cup,
rigid, double action




		Insert only	Complete instrument
340mm	5 Charr.	351-116-535	351-116-335
	7 Charr.	351-116-537	351-116-337
400mm	5 Charr.	351-116-505	351-116-305
	7 Charr.	351-116-507	351-116-307




Tenaculum Forceps, 1x1 teeth,
semi-rigid, double action

		Insert only	Complete instrument
340mm	5 Charr.	351-117-535	351-117-335
400mm		351-117-505	351-117-305



SHARK Tenaculum Grasping Forceps, semi-rigid, double action		Insert only	Complete instrument
340mm	5 Charr.	351-115-535	351-115-335
400mm		351-115-505	351-115-305



ACCARDI Harpoon Grasping Forceps, semi-rigid, double action		Insert only	Complete instrument
340mm	5 Charr.	351-119-535	351-119-335
400mm		351-119-505	351-119-305



Biopsy Forceps, rigid, double action		Insert only	Complete instrument
340mm	5 Charr.	351-120-535	351-120-335
	7 Charr.	351-120-537	351-120-337
400mm	5 Charr.	351-120-505	351-120-305
	7 Charr.	351-120-507	351-120-307



Biopsy Punch Forceps, semi-rigid, single action		Insert only	Complete instrument
340mm	5 Charr.	351-126-535	351-126-335
400mm		351-126-505	351-126-305

SCISSORS


- Standard handle
- White color coded



351-100-500
Standard handle
for scissors
(recommended handle)




Scissors, sharp tips, single action		Insert only	Complete instrument
340mm	5 Charr.	351-118-535	351-118-335
	7 Charr.	351-118-537	351-118-337
400mm	5 Charr.	351-118-505	351-118-305
	7 Charr.	351-118-507	351-118-307



Scissors, blunt tips, rigid, single action		Insert only	Complete instrument
340mm	5 Charr.	351-100-535	351-100-335
	7 Charr.	351-100-537	351-100-337
400mm	5 Charr.	351-100-505	351-100-305
	7 Charr.	351-100-507	351-100-307


ASSEMBLY



STEP 1

Open the handle.
The ball socket pushes out.


Keep the handle open.



STEP 2

1 Opened jaw


2 Close the jaw - the ball head pushes out. Keep the jaw closed.



STEP 3

Push back the locking nut.


With both opened handle and jaws, connect the ball head and ball socket.



STEP 4

Insert the ball head into the socket.

Keep the handle opened.



STEP 5

Close the handle. Shaft and insert will be retracted.
Keep the handle closed.

Screw the locking nut onto the handle.



DONE

Check the connection for being tight.

The knob rotates 360° and the handle as well as the jaw opens and closes.



ELECTRODES

360mm working length



Bipolar Needle Electrode,
straight tip
351-360-005



Bipolar Needle Electrode,
90° angled tip
351-360-015




Bipolar Ball Electrode,
360mm working length
351-370-105



ELECTRODES

360mm working length



Monopolar Needle Electrode
351-380-360





Monopolar Ball Electrode
351-385-360



SNARE

340mm working length



Polypectomy Snare
253-570-005

HIGH FREQUENCY CABLES

4.5m length




253-610-001* with 5 mm plug for Erbe HF-unit

253-610-002* with 4 mm plug for Martin HF-unit

253-610-003* with 8 mm plug for Valleylab HF-unit

253-610-004* with 4 mm plug for older HF-units

700-501-453* SDS plug




253-400-601 for Erbe HF-unit

253-400-602 for Martin HF-unit


700-500-453 SDS plug

253-400-606 for BOWA








RZ HYBRID RESECTOSCOPY SYSTEM





Complete Instrument consisting of **1-5**

► only cable needs to be changed

- 1** RZ Endoscope
- 2** HYBRID Electrode
- 3** 24 Charr. Continuous Flow Sheath
- 4** Titanium Working Element
- 5** HF Cable



5




ONE set
TWO techniques


18.5
CHARR.

HYBRID

GRAZIA RESECTOSCOPY SYSTEM



Endoscopes, suitable for HYBRID 18.5 Charr. System



Ø
2.9 mm



0° DIRECTION OF VIEW

300mm working length
356mm overall length

351-829-000



12° DIRECTION OF VIEW

300mm working length
356mm overall length

351-829-012



30° DIRECTION OF VIEW

300mm working length
356mm overall length

351-829-030

Bipolar Working Elements for GYNECARE-VERSAPOINT Bipolar Electrodes



351-000-400
Passive Titanium Working Element for use with
RZ Resectoscope Sheaths 253-000-256/253-000-356
and 4mm, 12°/30° RZ Hysteroscopes
351-804-012/351-804-030



351-000-401
Active Titanium Working Element for use with
RZ Resectoscope Sheaths 253-000-256/253-000-356
and 4mm, 12°/30° RZ Hysteroscopes
351-804-012/351-804-030

OLYMPUS compatible





351-000-500
Passive Titanium Working Element for use with
OLYMPUS Resectoscope Sheaths and Telescopes



351-000-501
Active Titanium Working Element for use with
OLYMPUS Resectoscope Sheaths and Telescopes

WIRE LOOP ELECTRODES

- Autoclavable handle with monopolar connector
- Loop which can be ordered as a spare part.
- To be used with a 5.5 mm trocar
- Salvage of the cervix, no vaginal descensus.
- Faster recovery of the patient.



350-754-000 Wire Loop Electrode with Ø 100mm wire loop electrode

Spare Wire Loops (single use, non-sterile)


- 350-754-001 Ø 100mm
- 350-754-002 Ø 125mm
- 350-754-003 Ø 65mm

MICRO NEEDLE ELECTRODES

- Extremely delicate needles
- Can be retracted during surgery



Pulled out needle,
cutting function



Retracted needle,
safe handling



Micro Needle Electrode for monopolar cutting

350-752-320 320mm

350-752-420 420mm

COMBINED MICRO NEEDLE ELECTRODE

- for cutting and coagulation



Pulled out needle,
cutting & coagulation



Retracted needle,
coagulation only




Combined Micro Needle Electrode for cutting and coagulation

350-753-320 320mm

350-753-420 420mm

350-753-001 Spare Micro Needle (Quantity: 12)

HIGH TEMPERATURE LED/XENON LIGHT CABLES



COMPATIBILITY NOTE:


A fiber optic cable with a large diameter should not be connected to a smaller diameter scope.
For example: $\varnothing$ 4.8mm fiber optic cables connected to a $\varnothing$ 2.9mm may risk high temperature at the light entry of the endoscope and may risk patient safety.

360° ROTATABLE INSTRUMENTS

Dedicated Instruments for Urology and Hysteroscopy

The unique design and features of the rotating color coded instruments have been designed in cooperation with Dr Accardi. Multiple years of experience in using dedicated instruments have led to this design and functionality.


The automatic holding mechanism combined with the 360° rotation of the jaw and the light desgin will improve your comfort.




360° ROTATABLE INSTRUMENTS

Instruments for Urology

Flexible and rigid instruments of 600mm length and different designs with various jaws for many indications in Urology.




Knob Color		5 Charr.	4 Charr.	
Teal	rigid	253-950-005	253-950-004	Stone Grasping Forceps, double action jaws
	flexible	253-951-005	253-951-004	
Blue	rigid	253-952-005	253-952-004	Stone Grasping Forceps, large stones, double action jaws
	flexible	253-953-005	253-953-004	
Yellow	rigid	253-954-005	253-954-004	HAMADA Long Stone Grasping Forceps, double action jaws
	flexible	253-955-005	253-955-004	
Black	rigid	253-956-005	253-956-004	Biopsy Forceps, double action jaws
	flexible	253-957-005	253-957-004	
Red	rigid	253-958-005	253-958-004	Scissors, single action jaws
	flexible	253-959-005	253-959-004	



360° ROTATABLE INSTRUMENTS

Instruments for Hysteroscopy


Semi-rigid instruments with dedicated jaws for your surgical hysteroscopy. The color codes will help always having the perfect instruments for your clinical needs.




	Knob Color	400mm	340mm	
		351-116-305	351-116-335	Grasping Forceps with cup, serrated, double action
		351-118-305	351-118-335	Scissors, sharp tips, single action
		351-100-305	351-100-335	Scissors, blunt tips, single action
		351-126-305	351-126-335	Biopsy Punch Forceps, single action
		351-120-305	351-120-335	Biopsy Forceps, double action
		351-117-305	351-117-335	Tenaculum Forceps, 1x1 teeth, double action
		351-115-305	351-115-335	SHARK Tenaculum Grasping Forceps
		351-119-305	351-119-335	Accardi HARPOON Grasping Forceps, with barbs



OBSTETRICS & GYNECOLOGY



Scissors
Clamps and Forceps
Hysterectomy Forceps
Uterus Forceps
Tenaculum Forceps
Uterine / Placenta Forceps
Obstetrical Forceps
Myoma Screws
Curettes
Amniotome
Uterine Scoops/Curettes
Biopsy Forceps
Probes/Dilators
Retractors/Specula
Pelvimeter



MAYO NOBLE

Scissors, curved

170mm, 6²/₃" 120-049-170

190mm, 7¹/₂" 120-049-190



WALDMANN

Episiotomy Scissors, curved, serrated

180mm, 7" 120-231-180



BRAUN STADLER

Episiotomy Scissors

145mm, 5³/₄" 120-225-145

200mm, 8" 120-225-200



SIMS

Gynecological Scissors, straight, blunt/blunt

200mm, 8"	120-062-200
-----------	-------------

230mm, 9"	120-062-230
-----------	-------------

curved

200mm, 8"	120-063-200
-----------	-------------

230mm, 9"	120-063-230
-----------	-------------



SIEBOLD

Gynecological Scissors, s-curved

240mm, 9 1/2"	120-061-240
---------------	-------------



WERTHEIM

Preparation Scissors, curved

195mm, 7 3/4"	120-053-195
---------------	-------------

225mm, 8 7/8"	120-053-225
---------------	-------------



Parametrium Scissors, TC
230mm, 9"

A straight	120-537-001
B slightly curved	120-537-002
C strongly curved	120-537-003
D angled	120-537-004

265mm, 10^{1/2}"





A straight	120-536-001
B slightly curved	120-536-002
C strongly curved	120-536-003
D angled	120-536-004

Parametrium Scissors, TC
305mm, 12"

A straight	120-538-001
B slightly curved	120-538-002
C strongly curved	120-538-003
D angled	120-538-004

Parametrium Scissors, TC
350mm, 14"

A straight	120-539-001
B slightly curved	120-539-002
C strongly curved	120-539-003
D angled	120-539-004





MODEL USA

Umbilical Cord Scissors

105mm, 4¹/₃"

120-233-105



BUSCH

Umbilical Cord Scissors, lateral curved

160mm, 6¹/₄"

120-229-160




SCHUMACHER

Umbilical Cord Scissors

160mm, 6¹/₄"

120-227-160




FAURE

Peritoneal Forceps, 200mm, 8"

Fig. 1 140-551-200
slightly curved

Fig. 2 140-553-200
strongly curved



MIKULICZ

Peritoneal Forceps, strongly curved, 1x2 teeth

180mm, 7" 140-549-185

200mm, 8" 140-549-200

laterally curved, 1x2 teeth

180mm, 7" 140-549-180




MIKULICZ

Peritoneal Forceps, curved, 1x2 teeth

180mm, 7" 140-548-180

200mm, 8" 140-548-200

245mm, 9²/₃" 140-548-245




LONG

Hysterectomy Forceps, 1x2 teeth

190mm, 7¹/₂"

140-652-190




GARLAND

Hysterectomy Forceps, 1x2 teeth

210mm, 8¹/₄"

140-639-210



PHANEUF

Peritoneal Forceps, 1x2 teeth,
straight


215mm, 8¹/₂"

140-520-210

curved

210mm, 8¹/₄"

140-521-210




HEANEY

Hysterectomy Forceps, delicate, single tooth, curved

210mm, 8¹/₄"

140-612-210



HEANEY


Hysterectomy Forceps, single tooth, curved

210mm, 8¹/₄"

140-613-210

230mm, 9"

140-613-230



HEANEY


Hysterectomy Forceps, double tooth, curved

210mm, 8¹/₄"

140-615-210

230mm, 9"

140-615-230



BALLANTINE


Hysterectomy Forceps, straight, single tooth

210mm, 8¹/₄" 140-617-210

300mm, 11⁷/₈" 140-617-300

curved


210mm, 8¹/₄" 140-618-210



HEANEY KANTER

Hysterectomy Forceps, curved, single tooth


210mm, 8¹/₄" 140-619-210



HEANEY REZEK

Hysterectomy Forceps, straight, single tooth



210mm, 8¹/₄" 140-627-210



GUPTA FRANK

Forceps, s-curved, 80mm jaw length



240mm, 9 ¹ / ₂ "	140-707-240
270mm, 10 ² / ₃ "	140-707-270
320mm, 12 ¹ / ₂ "	140-707-320



WERTHEIM

Hysterectomy Forceps, serrated, curved





250mm, 9 ⁷ / ₈ "	140-697-250
--	-------------



WERTHEIM CULLEN

Hysterectomy Forceps, 90° angled

215mm, 8 ¹ / ₂ "	140-641-215
--	-------------



A Parametrium Forceps, straight

210mm, 8 ¹ / ₄ "	140-700-220
250mm, 9 ⁷ / ₈ "	140-700-250
305mm, 12"	140-700-305
355mm, 14"	140-700-350

B slightly curved

210mm, 8 ¹ / ₄ "	140-701-215
240mm, 9 ¹ / ₂ "	140-701-245
305mm, 12"	140-701-305
355mm, 14"	140-701-350

C strongly curved

210mm, 8 ¹ / ₄ "	140-702-210
240mm, 9 ¹ / ₂ "	140-702-245
305mm, 12"	140-702-305
355mm, 14"	140-702-350

D angled

210mm, 8 ¹ / ₄ "	140-703-205
240mm, 9 ¹ / ₂ "	140-703-245
305mm, 12"	140-703-305
355mm, 14"	140-703-350



A straight, angled shanks

210mm, 8 ¹ / ₄ "	141-700-220
240mm, 9 ¹ / ₂ "	141-700-250

B slightly curved, angled shanks


210mm, 8 ¹ / ₄ "	141-701-215
240mm, 9 ¹ / ₂ "	141-701-245

C strongly curved, angled shanks

210mm, 8 ¹ / ₄ "	141-702-210
240mm, 9 ¹ / ₂ "	141-702-245

D angled, angled shanks

210mm, 8 ¹ / ₄ "	141-703-210
245mm, 9 ² / ₃ "	141-703-245




NOTO

Uterus Polypus Forceps

270mm, 10²/₃"

350-492-270




KELLY

Uterus Polypus Forceps

320mm, 12¹/₂"

350-490-320




GROSS MAIER

Swab Forceps, curved

220mm, 8²/₃"

150-031-220




BARRETT

Tenaculum Forceps

180mm, 7"

350-426-180




STAUDE MOORE

Tenaculum Forceps

170mm, 6²/₃"

350-422-170





POZZI

Tenaculum Forceps

260mm, 10¹/₄"

350-420-260




ADAIR

Tenaculum Forceps

190mm, 7 1/2"

350-426-190




JARCHO

Tenaculum Forceps

200mm, 8"

350-427-200




CZERNY

Tenaculum Forceps, straight

205mm, 8 1/8"

350-480-205



CORSON

Grasping Forceps, 330mm, 13"

7mm jaw width


350-755-007

11mm jaw width

350-755-011

17mm jaw width

350-755-017




SCHROEDER

Tenaculum Forceps

250mm, 9⁷/₈"

350-400-250




BRAUN

Tenaculum Forceps

250mm, 9⁷/₈"

350-410-250




DUPLAY

Tenaculum Forceps

280mm, 11"

350-428-280




HULKA

Uterine Manipulator Tenaculum Forceps

290mm, 11 1/2"

350-302-290



SCHROEDER


Tenaculum Forceps, 2 teeth, 240mm, 9 1/2"

straight

350-440-240

curved

350-441-240





TEALE

Tenaculum Forceps, laterally curved, 3x4 teeth

230mm, 9"

350-443-230





IOWA

Vulsellum Forceps

260mm, 10¹/₄"

350-445-260



JACOBS


Tenaculum Forceps, 2x2 teeth, 210mm, 8¹/₄"

straight

350-470-210

curved

350-471-210




GREEN ARMYTAGE

Caesarean Forceps, straight,
15mm wide serrated jaw

210mm, 8¹/₄"

340-150-210



MUSEUX

Forceps, 2x2 teeth, 8mm jaw
240mm, 9¹/₂" working length

straight

350-452-240

curved

350-453-240

10mm jaw

straight

350-454-240

curved

350-455-240



DARTIGUES

Uterine Elevating Forceps

270mm, 10²/₃"

350-493-270



SOMER

Uterine Elevating Forceps

230mm, 9"

350-493-230



COLLIN

Uterine Forceps

250mm, 9⁷/₈"

350-493-250



SAENGER

Placenta and Ovum Forceps,
270mm, 10²/₃"

straight	340-798-270
curved	340-799-270



Fig. 1




Fig. 2



WINTER

Placenta and Ovum Forceps, Fig. 2

280mm, 11"	340-799-280
------------	-------------

Placenta and Ovum Forceps, 280mm, 11"

Fig. 1	340-800-280
--------	-------------

Fig. 2	340-802-280
--------	-------------

Fig. 1 curved	
280mm, 11"	340-801-280



BIERER

Ovum Forceps, curved, 16mm jaw,
280mm, 11"

without ratchet	350-496-330
with ratchet	350-497-330




WRIGLEY

Obstetrical Forceps

280mm, 11"

340-011-280



SIMPSON

Obstetrical Forceps

300mm, 11^{7/8}"

340-013-300



SIMPSON BRAUN

Obstetrical Forceps

330mm, 13"

340-012-330

360mm, 14"

340-012-360

**SIMPSON LUIKART**

Obstetrical Forceps

360mm, 14"

340-029-360

**DE LEE**

Obstetrical Forceps

300mm, 11^{7/8}"

340-018-300

360mm, 14"

340-018-360



ELLIOT

Obstetrical Forceps

380mm, 15"

340-025-380




KIELLAND LUIKART

Obstetrical Forceps

400mm, 15^{3/4}"

340-016-400




KIELLAND

Obstetrical Forceps

410mm, 16^{1/8}"

340-014-410



DOYEN

Myoma screw

170mm

working length

350-757-190




DOYEN

Myoma screw

170mm working
length

350-759-170



Simple IUD hook with ring handle, malleable
250mm, 9^{7/8}" 350-766-001






Double IUD Extractor
270mm, 10^{2/3}" 350-766-010



SAUNDERS
IUD Removing Hook
250mm, 9^{7/8}" 350-766-100



EMMET
Fistula Hook, 220mm, 8^{2/3}"
Fig. 1 350-765-001
Fig. 2 350-765-002
Fig. 3 350-765-003
Fig. 4 350-765-004
Fig. 5 350-765-005






NOVAK

Biopsy Curette, 250mm, 9 7/8"

Ø 1mm	350-601-230
Ø 2mm	350-602-230
Ø 3mm	350-603-230
Ø 4mm	350-604-230




A Aspirating Curette, 290mm, 11 1/2"

3.0mm	350-650-030
4.5mm	350-650-045
6.0mm	350-650-060
8.0mm	350-650-080
10.0mm	350-650-100
12.0mm	350-650-120
14.0mm	350-650-140

B Coupling Piece

350-650-000




RANDALL

Suction Biopsy Curette, 240mm, 9 1/2"

Ø 4mm

350-602-240



KEYORKIAN

Suction Biopsy Curette, Ø 4mm,
without basket, rigid

300mm, 11 7/8"

350-599-300

with basket, rigid


300mm, 11 7/8"

350-600-300

with basket, malleable

300mm, 11 7/8"

350-601-300



MOD. MÜNCHEN

Suction Biopsy Curette, Ø 4mm

270mm, 10 2/3"

350-605-270



BEACHAM

Amniotome

260mm, 10¹/₄"

350-429-260




BRAUN

Uterine Depressor

280mm, 11"




350-399-280



SIMON

Uterine Scoop, sharp, 240mm, 9¹/₂"



Fig. 1, 7.0mm	350-610-001
Fig. 2, 8.5mm	350-610-002
Fig. 3, 9.5mm	350-610-003
Fig. 4, 11.5mm	350-610-004
Fig. 5, 13.0mm	350-610-005
Fig. 6, 14.5mm	350-610-006



SCHROEDER

Uterine Scoop, sharp, 300mm, 11⁷/₈"

Fig. 1, 6.8mm	350-612-001
Fig. 2, 8.5mm	350-612-002
Fig. 3, 10.0mm	350-612-003
Fig. 4, 11.5mm	350-612-004
Fig. 5, 13.0mm	350-612-005
Fig. 6, 14.5mm	350-612-006



GOURDET

Uterine Scoop, blunt, 280mm, 11"

Fig. 1, 6mm	350-613-006
Fig. 2, 8mm	350-613-008
Fig. 3, 10mm	350-613-010
Fig. 4, 12mm	350-613-012
Fig. 5, 14mm	350-613-014



THOMAS

Uterus Curette, malleable, blunt, 280mm, 11"

Fig. 1, Ø 10mm	350-919-001
Fig. 2, Ø 11mm	350-919-002
Fig. 3, Ø 12mm	350-919-003
Fig. 4, Ø 13mm	350-919-004
Fig. 5, Ø 14mm	350-919-005
Fig. 6, Ø 15mm	350-919-006



SIMS

Uterus Curette, 285mm, 11^{1/4}"

rigid, sharp

rigid, blunt

malleable, sharp

malleable, blunt



Fig. 00, 5mm



Fig. 0, 6mm



Fig. 1, 7mm



Fig. 2, 8mm

350-630-005

350-630-006

350-630-007

350-630-008

350-926-050

350-926-060

350-926-070

350-926-080

350-631-005

350-631-006

350-631-007

350-631-008

350-927-050

350-927-060

350-927-070

350-927-080



Fig. 3, 9mm



Fig. 4, 11mm



Fig. 5, 12mm




Fig. 6, 14mm

rigid, sharp

350-630-009

350-630-011

350-630-012

350-630-014

rigid, blunt

350-926-090

350-926-110

350-926-120

350-926-140

malleable, sharp

350-631-009

350-631-011

350-631-012

350-631-014

malleable, blunt

350-927-090

350-927-110

350-927-120

350-927-140




Fig. 7, 15mm



Fig. 8, 16mm



Fig. 9, 19mm




Fig. 10, 20mm

rigid, sharp

350-630-015

350-630-016

350-630-019

350-630-020

rigid, blunt

350-926-150

350-926-160

350-926-190

350-926-200

malleable, sharp

350-631-015

350-631-016

350-631-019

350-631-020

malleable, blunt

350-927-150

350-927-160

350-927-190

350-927-200



RECAMIER

Uterus Curette, 310mm, 12^{1/4}"

	Fig. 00, 5mm	Fig. 0, 6mm	Fig. 1, 7mm	Fig. 2, 8mm
rigid, sharp	350-620-005	350-620-006	350-620-007	350-620-008
rigid, blunt	350-619-005	350-619-006	350-619-007	350-619-008
malleable, sharp	350-621-005	350-621-006	350-621-007	350-621-008
malleable, blunt	350-622-005	350-622-006	350-622-007	350-622-008



Fig. 00, 5mm



Fig. 0, 6mm



Fig. 1, 7mm



Fig. 2, 8mm



Fig. 3, 9mm



Fig. 4, 11mm



Fig. 5, 12mm




Fig. 6, 14mm

rigid, sharp	350-620-009	350-620-011	350-620-012	350-620-014
rigid, blunt	350-619-009	350-619-011	350-619-012	350-619-014
malleable, sharp	350-621-009	350-621-011	350-621-012	350-621-014
malleable, blunt	350-622-009	350-622-011	350-622-012	350-622-014




Fig. 7, 15mm



Fig. 8, 16mm




Fig. 9, 19mm





Fig. 10, 20mm

rigid, sharp	350-620-015	350-620-016	350-620-019	350-620-020
rigid, blunt	350-619-015	350-619-016	350-619-019	350-619-020
malleable, sharp	350-621-015	350-621-016	350-621-019	350-621-020
malleable, blunt	350-622-015	350-622-016	350-622-019	350-622-020



BUMM

Uterus Curette, 300mm, 11^{7/8}"

Fig. 1, 30mm

Fig. 2, 35mm

Fig. 3, 40mm

rigid, sharp

350-920-300

350-920-350

350-920-400

rigid, blunt

350-926-300

350-926-350

350-926-400

malleable, sharp

350-928-300

350-928-350



350-928-400

malleable, blunt

350-929-300

350-929-350

350-929-400




BLAKE

Biopsy Curette, sharp/blunt

270mm, 10^{2/3}"

350-614-270




WALLICH

Placenta / Ovum Scoop

420mm, 16^{1/2}"

340-510-420



PESTALOZZA

Uterus Curette

300mm, 11^{7/8}"

340-500-300



GREENE

Uterus Curette, 280mm, 11"

10mm

350-365-010

12mm

350-365-012

14mm

350-365-014



MILAN

Curette

290mm, 11^{1/2}"

350-605-007



FAURE

Biopsy Forceps

210mm, 8¹/₄"

350-843-210



SCHUBERT

Biopsy Forceps

210mm, 8¹/₄"

350-844-210

290mm, 11¹/₂"

350-844-290




SCHUBERT

Biopsy Forceps

260mm, 10¹/₄"

350-847-260



THOMS GAYLOR

Biopsy Forceps, curved

220mm, 8 ² / ₃ "	350-857-220
--	-------------


240mm, 9 ¹ / ₂ "	350-857-240
--	-------------



VAN DOREN

Biopsy Forceps


267mm, 10 ¹ / ₂ "	350-855-267
---	-------------



ALEXANDER

Biopsy Forceps

280mm, 11"	310-735-280
------------	-------------




KERVORKIAN

Uterine Biopsy Forceps, 9x3mm

240mm, 9^{1/2}"

310-730-240



TISCHLER MORGAN


Biopsy Forceps, 7x3mm

200mm, 8"

310-733-001

250mm, 9^{7/8}"

310-735-001



MINI TOWNSEND


Biopsy Forceps, 4.2 x 2.3mm jaw

200mm, 8"

310-733-002

250mm, 9^{7/8}"

310-735-002




BABY TISCHLER

Biopsy Forceps, 4.2 x 2.3mm jaw

200mm, 8" 310-733-003


250mm, 9⁷/₈" 310-735-003



TISCHLER KEVORKIAN

Biopsy Forceps, 8 x 3mm jaw

250mm, 9⁷/₈" 310-731-250



KEVORKIAN PACIFIC


Biopsy Forceps, 9 x 3mm jaw

200mm, 8" 310-732-200

250mm, 9⁷/₈" 310-732-250

8 x 3mm jaw

200mm, 8" 310-731-200



Detachable Biopsy Forceps "Roto Fit"



UNI TISCHLER

Insert only, 3 x 7mm, 230mm, 9"

A straight 310-800-100

B angled up 310-801-100

C angled down 310-802-100

complete instrument with handle

A straight 310-800-000

B angled up 310-801-000

C angled down 310-802-000

UNI TOWNSEND

Insert only, 2.3 x 7mm, 230mm, 9"

D straight 310-803-100

E angled up 310-804-100

F angled down 310-805-100

complete instrument with handle

D straight 310-803-000

E angled up 310-804-000

F angled down 310-805-000

UNI KERVOKIAN

Insert only, 3 x 9mm, 230mm, 9"

G straight 310-806-100

complete instrument with handle

G straight 310-806-000

UNI OVAL

Insert only, 3 x 5.5mm, 230mm, 9"

H straight 310-807-100

complete instrument with handle

H straight 310-807-000

I Handle only


310-800-600



BURKE

Biopsy Forceps, 5 x 3mm jaw


200mm, 8"	310-732-000
250mm, 9 ⁷ / ₈ "	310-735-004



EPPENDORFER

Biopsy Forceps, 4x8mm

200mm, 8"	310-760-200
230mm, 9"	310-760-230




SCHUMACHER

Biopsy Forceps,
Fig. 1

240mm, 9 ¹ / ₂ "	310-729-001
--	-------------

Fig. 2

240mm, 9 ¹ / ₂ "	310-729-002
--	-------------



Transobturator Sling needle
pair of right and left curved

260mm, 10¹/₄" 350-101-260



MIYAZAKI

Ligature Carrier, large hook 0.75"

200mm, 8" shaft 180-160-000



MIYAZAKI

Speculum with notch

180-160-010



SIMS

Uterine Probe, graduated, Ø 4mm,
330mm, 13"

rigid 350-395-320

malleable 350-396-320



SIMPSON

Uterine Probe, malleable, Ø 4mm

290mm, 11¹/₂" 350-396-290



Uterine manipulator, complete set
(350-660-111 - 350-660-116)

350-660-100



Adaptor Ø 8 x 80 mm

350-660-111



Adaptor Ø 10 x 100 mm

350-660-112



Adaptor Ø 3 x 45 mm

350-660-113



Adaptor Ø 3 x 55 mm

350-660-114



Adaptor Ø 3.2 x 35 mm

350-660-115



Adaptor Ø 3.2 x 45 mm

350-660-116



COHEN

Uterus Cannula
complete with 2 interchangeable acorns

350-750-000

Acorns

small 350-750-001

large 350-750-002



HEGAR

Uterine Dilator, double-ended

1.0/2.0mm	350-382-020	15.0/16.0mm	350-382-160
3.0/4.0mm	350-382-040	17.0/18.0mm	350-382-180
5.0/6.0mm	350-382-060	19.0/20.0mm	350-382-200
7.0/8.0mm	350-382-080	21.0/22.0mm	350-382-220
9.0/10.0mm	350-382-100	23.0/24.0mm	350-382-240
11.0/12.0mm	350-382-120	25.0/26.0mm	350-382-260
13.0/14.0mm	350-382-140		

Set*, 8 sizes from 3/4 to 17/18 mm
in canvas roll

350-383-001



HEGAR

Uterine Dilator with slope handle

1.0mm	350-380-010	12.5mm	350-380-125	24.0mm	350-380-240
1.5mm	350-380-015	13.0mm	350-380-130	25.0mm	350-380-250
2.0mm	350-380-020	13.5mm	350-380-135	26.0mm	350-380-260
2.5mm	350-380-025	14.0mm	350-380-140	27.0mm	350-380-270
3.0mm	350-380-030	14.5mm	350-380-145	28.0mm	350-380-280
3.5mm	350-380-035	15.0mm	350-380-150	29.0mm	350-380-290
4.0mm	350-380-040	15.5mm	350-380-155	30.0mm	350-380-300
4.5mm	350-380-045	16.0mm	350-380-160		
5.0mm	350-380-050	16.5mm	350-380-165		
5.5mm	350-380-055	17.0mm	350-380-170		
6.0mm	350-380-060	17.5mm	350-380-175		
6.5mm	350-380-065	18.0mm	350-380-180		
7.0mm	350-380-070	18.5mm	350-380-185		
7.5mm	350-380-075	19.0mm	350-380-190		
8.0mm	350-380-080	19.5mm	350-380-195		
8.5mm	350-380-085	20.0mm	350-380-200		
9.0mm	350-380-090	21.0mm	350-380-210		
9.5mm	350-380-095	22.0mm	350-380-220		
10.0mm	350-380-100	23.0mm	350-380-230		
10.5mm	350-380-105				
11.0mm	350-380-110				
11.5mm	350-380-115				
12.0mm	350-380-120				

* illustration does not show the complete set



PRATT

Uterine Dilators, French Scale 10³/₈"

13 - 15	350-392-015
17 - 19	350-392-019
21 - 23	350-392-023
25 - 27	350-392-027
29 - 31	350-392-031
33 - 35	350-392-035
37 - 39	350-392-039
41 - 43	350-392-043

Set of 8 pieces in canvas roll*

350-392-000



HANK

Uterine Dilator


Fig. 9/10 4.5/5.0mm	350-390-010
Fig. 11/12 5.5/6.0mm	350-390-012
Fig. 13/14 6.5/7.0mm	350-390-014

Fig. 15/16 7.5/8.0mm	350-390-016
Fig. 17/18 8.5/9.0mm	350-390-018
Fig. 19/20 9.5/10.0mm	350-390-020

Set, 6 pieces in canvas roll*

350-391-002

* illustration does not show the complete set



JACKSON

Vaginal Speculum

75 x 38mm	350-206-001
90 x 38mm	350-206-002
100 x 38mm	350-206-003



EASTMAN

Vaginal Speculum

75 x 38mm	350-207-001
90 x 38mm	350-207-002
100 x 38mm	350-207-003



DOYEN

Vaginal Speculum, slightly concave

55 x 35mm	350-290-001
85 x 35mm	350-290-002
115 x 35mm	350-290-003
55 x 45mm	350-290-004
85 x 45mm	350-290-005
115 x 45mm	350-290-006
160 x 45mm	350-290-007
55 x 60mm	350-290-008
85 x 60mm	350-290-009
115 x 60mm	350-290-010
160 x 60mm	350-290-011



HEANEY

Retractor, 270mm, 10²/₃"

73 x 28mm	170-459-001
45 x 28mm	170-459-002
82 x 28mm	170-459-003



HEANEY SIMON

Retractor with horn handle


260mm, 10 ¹ / ₄ "	170-460-260
---	-------------



BREISKY

Vaginal Speculum

100 x 20mm	350-320-020	130 x 20mm	350-322-020
100 x 25mm	350-320-025	130 x 25mm	350-322-025
100 x 30mm	350-320-030	130 x 30mm	350-322-030
100 x 35mm	350-320-035	130 x 35mm	350-322-035
100 x 40mm	350-320-040	130 x 40mm	350-322-040
		155 x 40mm	350-324-040



SEIDL

Speculum for children
170mm, 6^{2/3}"

80 x 8mm	350-254-001
80 x 10mm	350-254-002
90 x 14mm	350-254-003



KRISTELLER

Vaginal Speculum, for children

70 x 15mm	350-210-000
for adults	
70 x 27mm	350-210-001
80 x 30mm	350-210-002
90 x 36mm	350-210-003



KRISTELLER

Retractor

Fig. 0	70 x 10mm	350-211-000
Fig. 1	75 x 23mm	350-211-001
Fig. 2	85 x 26mm	350-211-002
Fig. 3	100 x 32mm	350-211-003



KRISTELLER

Vaginal Speculum Set

Fig. 0	350-212-000
Fig. 1	350-212-001
Fig. 2	350-212-002
Fig. 3	350-212-003

**KRISTELLER**

Vaginal Speculum

Fig. 1	110 x 27mm	350-230-001
Fig. 2	110 x 30mm	350-230-002
Fig. 3	110 x 36mm	350-230-003

**KRISTELLER**

Retractor

Fig. 1	115 x 18mm	350-231-001
Fig. 2	115 x 22mm	350-231-002
Fig. 3	115 x 28mm	350-231-003

**KRISTELLER**

Vaginal Speculum Set

Fig. 1	350-232-001
Fig. 2	350-232-002
Fig. 3	350-232-003



SIMS

Vaginal Speculum, small

65 x 25mm	350-200-001
-----------	-------------

medium

70 x 30mm	350-200-002
-----------	-------------

large

85 x 35mm	350-200-003
-----------	-------------



RZ

Lateral Retractor

60 x 22mm	350-340-000
-----------	-------------



KOGAN

Endospeculum, screw ratchet, 240mm, 9 1/2"

5mm	350-361-001
-----	-------------

4mm	350-361-002
-----	-------------

3mm	350-361-003
-----	-------------



GRAVE

Vaginal Speculum

75 x 20mm	350-130-001
95 x 35mm	350-130-002
115 x 35	350-130-003



PEDERSON

Vaginal Speculum


75 x 13mm	350-140-001
100 x 22mm	350-140-002
120 x 25mm	350-140-003



CUSCO

Vaginal Speculum

75 x 32mm	350-100-001
85 x 35mm	350-100-002
100 x 35mm	350-100-003



CUSCO

Vaginal Speculum, Swiss pattern

80 x 22mm	350-102-001
100 x 25mm	350-102-002
105 x 27mm	350-102-003



SEMN

Vaginal Speculum, 100mm blade

17-20mm wide blade	350-120-001
25-30mm wide blade	350-120-002



COLLIN

Vaginal Speculum, detachable

65 x 20mm	350-110-001
85 x 30mm	350-110-002
100 x 35mm	350-110-003
100 x 40mm	350-110-004



AUVARD

Vaginal Speculum, fixed weight,
1.14kg, 45° angled

70 x 42mm	350-309-000
83 x 38mm	350-310-000
102 x 45mm	350-311-000



STEINER ANVARD

Vaginal Speculum, fixed weight, 1.14kg

150 x 34mm	350-311-140
------------	-------------



SCHERBACK

Vaginal Speculum, complete set
incl. weight 400g, handle, blades Fig. 1 - 5

350-300-000



Blades	
80x30mm	330-301-001
85x35mm	330-301-002
85x40mm	330-301-003
55x40mm	330-301-004
70x45mm	330-301-005



Weight	
400g	350-302-001
800g	350-302-002
Handle	
	350-303-000




GELPI

Perineal Retractor, pediatric size

90mm, 3 1/2" 170-802-090

110mm, 4 1/3" 170-802-110

135mm, 5 1/3" 170-802-135



GELPI

Perineal Retractor

180mm, 7" 170-802-180

with ball stops

190mm, 7 1/2" 170-803-180



RIGBY

Perineal Retractor, 21x70mm blades

180mm, 7" 170-806-180



O'SULLIVAN O'CONNOR

Vaginal Retractor Set

350-336-000




RIGBY

Perineal Retractor, 21x70mm blades

180mm, 7"

170-806-180





MARTIN

Pelvimeter, cm/inch scale

340mm, 13 1/2"

340-103-340





C2 CAMERA & LIGHT SOURCE COMBINATION

This integrated FullHD-Camera and LED Light Source is a combined 2-in-1 system. The concept is specially designed for the use in the outpatient market.

The C2 support rigid scopes and fiber-endoscopes and with the small footprint it provides an economical and flexible solution for endoscopic diagnostic and surgical interventions in Gynecology, Urology or ENT.

Features:

- Compact and integrated 2-in-1 design
- Versatile use for rigid scopes and fiberoptics
- LED lamp – exceptionally long lifespan
- Low noise ventilation system
- Easy storage and documentation of videos and pictures
- Optional TV couplers with fixed focal length or with zoom available




FULL HD



Integrated recording functions



Rigid and fiber endoscopes supported



Easy to use



250-015-030
Zoom TV Adapter

DELIVERY CONTAINS:

- Main Unit
- Camera head
- HDMI cable
- Mains cable
- Operating manual

C2 Camera and Light Source Combination
300-002-200

Technical data

Line voltage	100-240 VAC	
Power consumption	60VA	
Image sensor	1/2" MOS	
Resolution	1920 x 1080 (FullHD)	
Index ratio	16:9	
Lamp	LED (50W comparable)	
Video outputs	1x HDMI	1x USB (A)
	1x DVI	1x Remote IN (Footswitch)
	1x SDI	1x Ri45 (Service)
	1x VIDEO	
Dimensions	Control unit: 295 x 128 x 279mm (WxHxD)	
	Weight control unit: approx. 3kg	
	Weight camera head: 180g	

TV Coupler, Endoscope and Light Guide cable not included and can be ordered seperately.
NOTE: Endoscopes smaller than Ø 2.9mm require 250-015-030 Zoom TV Adapter.